

Aircare's wall-mounted drying rack provides safe and convenient storage for up to 24 test tubes. After tubes have been washed, they are placed on the pegs, drip out their water, air dry, and await reuse.

Made from vacuum-formed ABS with PVC pegs, the rack is easy to clean, has no sharp edges, and has an integral drainage basin and outlet. The standard rack size is 600mm x 600mm.

Specifications

Materials	<ul style="list-style-type: none"> ▶ Vacuum-formed 3mm ABS body ▶ 10mm diameter PVC pegs ▶ PVC drain connector
External Dimensions	600mm Wide x 600mm in Height

Options

- ▶ Overall height, overall width, and number of tube drying rods can be varied to suit.

